

Data Sheet:

Safety Sentry Pull-String Model 91230

Features

- Easy-to-use, low cost, pull-string monitor
- Single person monitoring device that alarms when unsupervised resident movement occurs
- Monitor attaches to a bed or chair; locking Safety Klip is attached to the resident
- Alarm sounds when magnetic connection is broken between magnet and monitor, indicating resident movement

Safety Klip

This device is locked onto an article of the resident's clothing. The unique plastic tab makes it difficult for the resident to detach without the assistance of a caregiver.



Shown with Optional
Cushion Protector
Model 0141-409



Universal Bed Bracket
Model 0707-578